

## Noseband and Chest Strap Cover Instructions

Thank you for your recent interest in making nose band and chest strap covers for SPANA. We are always grateful to receive more of these covers to help us with our ongoing work.

In the countries where SPANA works, we often see wounds on the noses and chests of horses, mules and donkeys, due to badly designed and ill-fitting harnessing. While in the long term we aim to help change the methods of harnessing used overseas we can also help to relieve a lot of the suffering by distributing noseband and chest strap covers kindly made for us by supporters like you.


### Instructions for making noseband and harness covers

- 1) Washable Nylon fur fabric preferred (but any durable soft material will do) the brighter and fluffier the better!
- 2) Noseband cover: size 11" x 4" Chest strap cover: size 24" x 7"
- 3) One strip of Velcro to be attached to the inside length; another strip to be attached to the outside/fluffy side of the length of the fabric. (this is better stitched than glued since the glue deteriorates easily in the heat)
- 4) The noseband can then be attached to the noseband or chest strap by rolling it around the existing noseband/chest strap and securing it with the Velcro.
- 5) Post any finished items to our SPANA head office, 14 John Street, London, WC1N 2EB

////////VELCRO////////

RIGHT SIDE OF FABRIC

(VELCRO ATTACHED TO UNDERSIDE)



These noseband and chest strap covers are easy for us to take overseas as they are light and can be packed around other essential medicines and equipment, and they are very popular with owners.

We constantly need them, since they do of course wear out.

SPANA is totally reliant on the generosity of its supporters to carry out its work in helping improve the lives of working animals and their owners around the world so thank you for considering helping SPANA in this way.

[www.spana.org](http://www.spana.org)